Privileges

The following table describes the user role privileges.

Privileges

PRIVILEGE NAME	DESCRIPTION
Change Lab Preferences	This privilege allows you to configure default units and names within projects and tests. Without this privilege, the application disables:
	• Preferences > Default Names
	• Preferences > Units > Units Management
Create Diagnostic File	This privilege allows the user to create a diagnostic file, which is used by MTS for diagnostic purposes. Without this privilege, the application disables: Tools > Create Diagnostic Package .
Edit Procedure	This privilege allows the user to define new tests and modify existing tests.
	Without this privilege, the application changes in several ways:
	 The Project Explorer does not show the Test Definition node under each test.
	• File > New > Test is disabled. Without this privilege, the only way to create a new test is from a template or existing project.
	• The application does not show the Toolbox tab.
	• The application does not show the Properties tab.
	• The Display menu does not show the Toolbox or Properties items.
Launch the Fatigue Analyzer Application	Start the Fatigue Analyzer application. Without this privilege, the application issues a message when you attempt to open the application.
Launch the Fracture Analyzer Application	Start the Fracture Analyzer application. Without this privilege, the application issues a message when you attempt to open the application.
Launch the Multipurpose Application	Start the Multipurpose application. Without this privilege, the application issues a message when you attempt to open the application.
Launch the Reporter Application	Start the Reporter application. Without this privilege, the application issues a message when you attempt to open the application.
User Management	View information related to users, roles, and privileges. Without this privilege, the application disables Preferences > Settings > User Management .
View Reports	View reports generated by an activity in a test procedure. Without this privilege, the application disables a right-click context menu View Report option.